

Luca Bloom

Instructions for Use, Safety Warnings & Terms of Sale — Australia

IMPORTANT — READ THIS ENTIRE DOCUMENT BEFORE USING THIS PRODUCT. IT CONTAINS SAFETY WARNINGS AND AFFECTS YOUR LEGAL RIGHTS.

We recommend saving or printing this document for future reference, including for emergency use.

By purchasing and/or using this product, you acknowledge that you have read, understood, and agree to be bound by all terms set forth below. If you do not agree, do not use this product.

Nothing in this document removes or reduces any right you have under the Australian Consumer Law.

1. Age Requirement

- **You must be at least 18 years of age to purchase or use this product. This product is not intended for, marketed to, or sold to anyone below that age.**
- The age of majority is 18 throughout Australia, in every state and territory.
- By purchasing this product, you represent and warrant that you meet this age requirement. Misrepresentation of your age is a material breach of these terms.

2. What This Product Is and Where It Is Used

Description

Luca Bloom is a hand-operated mechanical device made of POM (acetal) resin. It consists of a body, two slender arms, a trapezoidal lead screw with a turning knob, and a keyhole-shaped release plate.

The side of the body is engraved with the word **CLOSE** next to an arrow. This marking indicates the direction in which the knob must be turned to bring the two arms back together and reduce the expansion. It refers to the rotation of the knob only; it does not refer to the release plate, which is a separate part operated in a different way.

Step-by-step operating instructions are given separately in the illustrated instruction guide provided with this product. Read that guide together with this document before first use.

Application area — read this carefully

This product is intended for external use by adult males on the foreskin (prepuce) of the penis. The tips of the two arms are inserted a short distance into the opening of the foreskin. Turning the knob slowly separates the arms, which gently widens that opening for stretching purposes.

Throughout this document, “the application area” means the foreskin, the opening of the foreskin, the glans, and any adjacent tissue that the device contacts.

- **For external use only. Never insert the arms into the urethra (the urinary opening).**
- **Never insert the arms further than needed to make contact with the inner edge of the foreskin opening.**
- **This product is not intended for use on any other part of the body, or by anyone other than the individual purchaser.**

Not a therapeutic good

- This product is NOT a medical device or therapeutic good within the meaning of the Therapeutic Goods Act 1989 (Cth). It has not been assessed, approved, or evaluated by the Therapeutic Goods Administration (TGA), and it is not entered on the Australian Register of Therapeutic Goods.

- This product is not intended to diagnose, treat, cure, mitigate, or prevent any disease or medical condition. We make no medical or therapeutic claim, and no particular or permanent result is guaranteed.
- If you have any medical concern regarding the application area, consult a registered medical practitioner before using this product. This product is not a substitute for professional medical care.

3. Assumption of Risk

USE OF THIS PRODUCT INVOLVES INHERENT RISK OF SERIOUS PHYSICAL INJURY TO THE FORESKIN, GLANS, AND SURROUNDING TISSUE. By using this product, you knowingly and voluntarily assume the risks associated with its use, including risks arising from your own decisions regarding whether, when, how long, and how forcefully to use it.

Improper use, excessive force, over-extension, prolonged use, or unsanitary use may result in serious injury, including but not limited to: laceration or tearing of the foreskin, bruising, swelling, restricted or obstructed blood flow to the glans, nerve damage, loss of sensation, infection, adhesion, permanent scarring, and tissue necrosis. Some injuries may be permanent and may require surgical intervention.

You are responsible for reading, understanding, and following every instruction, contraindication, and warning in this document.

This section does not exclude, restrict or modify the consumer guarantees under the Australian Consumer Law, or any other liability that cannot be excluded under applicable law. See Sections 9 and 10.

4. Contraindications — DO NOT USE THIS PRODUCT IF ANY OF THE FOLLOWING APPLY

- You are unable to independently and safely judge, attach, operate, or remove the device.
- You have a severe metal allergy.
- You take anticoagulant (blood-thinning) medication, have any bleeding disorder, or have a severe vascular condition.
- You are undergoing treatment related to the application area, or you have reduced sensation due to a circulatory disorder, peripheral neuropathy, diabetes, or a similar condition.
- You have any wound, bleeding, inflammation, infection, discharge, rash, or other abnormality on or within the application area.
- **The foreskin is already adhered or fused to the glans. Forced expansion in this condition can cause severe tearing. Consult a medical practitioner instead.**
- You are under the influence of alcohol, or have taken pain medication, sleep aids, sedatives, or similar substances.
- Children or infants are nearby, or the product would be stored within their reach (choking and injury hazard).
- You intend to use it for more than 30 minutes in a single session, for more than one hour in total in any 24-hour period, or while sleeping. Each of these is prohibited.

5. Mandatory Safety Rules

Stop immediately and seek medical care

If you experience pain, congestion, swelling, numbness, discoloration, or any unexpected problem, stop using the product immediately and consult a medical practitioner promptly. In an emergency, call 000, or go to the nearest hospital emergency department.

Monitor the glans for colour change during use

To detect impaired blood flow, check the glans and foreskin regularly during use for discoloration (reddish-purple, darkening, or similar). Stop immediately if any abnormal discoloration appears.

Relieving pressure — the CLOSE marking

If you feel discomfort, pain, tightness, numbness, or any other unexpected sensation, the first thing to do is turn the knob in the direction of the arrow marked CLOSE on the side of the body. This brings the arms back together and reduces the pressure on the tissue.

Do this before attempting anything else. Do not pull the device, and do not try to force it off. In most cases, releasing the expansion in this way is enough to remove the device normally.

If the knob will not turn at all, do not force it and do not try to relieve pressure this way. Go directly to the release plate described below.

The release plate is also the correct step if the device still will not come off after the expansion has been released.

Emergency release — how to remove the device if it will not come off

DO NOT PULL THE DEVICE FORCEFULLY. Forcing it can cause serious injury.

This product includes a release plate — a flat, keyhole-shaped part. In its locked position, the narrow slot of the keyhole engages the shaft and holds the device in place. To release the device, push the release plate downward. This aligns the opening with the shaft so that the shaft can pass through it.

If the device still cannot be removed after operating the release plate, do not continue to force it. Seek medical attention promptly.

Never force the foreskin down

Never forcibly retract the foreskin when there is insufficient slack. Doing so can severely obstruct blood flow to the glans and cause serious tissue damage, including necrosis. If the retracted foreskin becomes trapped behind the glans and cannot be returned to its normal position, seek emergency medical care immediately.

Adhesion warning

If the inner surface of the foreskin is already adhered to the glans, forced expansion can cause severe lacerations. Do not attempt to force it open. Consult a medical practitioner.

No excessive expansion

Never expand to a degree that causes significant pain. Stop well before the point of pain.

Do not overload the arms

The arms are long and slender by design. Excessive force or load may break them, creating a risk of injury from fragments — including inside the application area. Do not apply excessive force during operation or cleaning.

If the device breaks during use

Stop immediately and remove the device. Recover all fragments. If a fragment remains in the application area, or if you cannot confirm that every fragment has been recovered, do not attempt to remove it yourself. Seek medical attention promptly.

Do not use a product that has broken, even if it appears to still function.

Lubricants, creams, and other products

After use, wash any residue off the device, particularly from the lead screw, so that the mechanism does not accumulate build-up. Avoid prolonged contact between the device and strongly acidic preparations, such as exfoliating products containing glycolic or salicylic acid.

Inspect before every use

Before every use, visually inspect the product — particularly the arms and the release plate — for cracks or damage. Do not use a damaged product. Take care not to pinch skin or hair during operation.

6. Material & Allergy Notice

- This product is made of POM (polyoxymethylene / acetal resin). All parts are POM.
- Because manufacturing takes place in a facility that also processes metal components, trace amounts of metal residue may be present on the product. Discontinue use immediately if you experience redness, irritation, or any other skin reaction.

Material composition note

The material manufacturer's published safety data for POM (acetal) resin identifies formaldehyde in connection with this resin, which may be present in trace amounts. This note is provided for your information. If you experience any reaction that you believe may be related, discontinue use and consult a medical practitioner.

7. Care, Storage & Disposal

Included items

Luca Bloom unit ×1; information card with QR code linking to the illustrated instruction guide and to this document ×1. The card is shared across markets; select your country at the address on the card.

Origin and finishing

Components are manufactured in China. The product is assembled, ultrasonically cleaned, and inspected in Japan before dispatch.

Before first use (required)

Machining oil and fine particulate residue from manufacturing may remain on the product. Wash thoroughly with a mild detergent and water, then dry completely before first use.

Routine care

Wash with a mild detergent after every use and dry completely before storage. If you disinfect the product, use isopropyl alcohol. Do not boil the product or expose it to hot water, as this may cause deformation.

Storage

Store at room temperature in a clean location away from direct sunlight, heat, and humidity. Keep out of reach of children.

Disposal

Dispose of the product in accordance with the plastic-waste regulations of your local council.

No modification

Do not modify, alter, or disassemble the product. Doing so voids the warranty described in Section 9 and constitutes misuse.

8. Resale & Transfer Prohibited

- For hygiene and safety reasons, resale, lending, gifting, or transfer of this product to any third party is strictly prohibited. This product is intended for use by the original purchaser only.
- Luca Bloom accepts no responsibility for any product obtained through an unauthorised channel, and the warranty in Section 9 does not apply to such products.

9. Australian Consumer Law & Warranty Against Defects

Your rights under the Australian Consumer Law

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

These consumer guarantees apply automatically and cannot be excluded, restricted, or modified by anything in this document. Where any provision of this document conflicts with them, the Australian Consumer Law prevails and that provision does not apply to you.

Warranty against defects

The benefits given to you by this warranty are in addition to other rights and remedies you have under a law in relation to the goods to which this warranty relates.

What this warranty covers: Luca Bloom warrants that the product will be free from manufacturing defects that render it inoperable upon arrival.

Warranty period: seven (7) days from the date of delivery.

How to claim: contact us by email or telephone using the details below within seven (7) days of delivery. State your order number and describe the defect. We will respond with instructions.

Remedy: at our option, a replacement product or a refund of the purchase price.

Who bears the expenses of claiming: Luca Bloom bears all expenses of claiming under this warranty, including any return postage. You will not be asked to pay any cost in order to make a claim.

Who gives this warranty

Luca Bloom

S-Building 3F, 2-1-19 Roppongi, Minato-ku, Tokyo, Japan

Telephone: +81 70-3842-6430

Telephone enquiries are answered in Japanese only, on weekdays from 10:00 to 17:00 Japan Standard Time. For enquiries in English, please use email or the platform messaging system; we answer in English by those channels.

Email: lucabloom65@gmail.com

What this warranty does not cover

Damage from dropping, misuse, improper operation, modification, disassembly, or normal wear; products obtained through unauthorised channels; or any personal injury. Personal injury is dealt with in Section 10 and is not subject to this warranty.

Returns

Because this is a personal hygiene product that comes into direct contact with the body, returns and exchanges for change of mind or personal preference cannot be accepted, except where a right of return is granted to you by law or by the platform through which you purchased.

10. Limitation of Liability

Consumer guarantees are not excluded

NOTHING IN THIS DOCUMENT EXCLUDES, RESTRICTS OR MODIFIES THE CONSUMER GUARANTEES UNDER THE AUSTRALIAN CONSUMER LAW, OR ANY OTHER LIABILITY THAT CANNOT BE EXCLUDED OR LIMITED UNDER APPLICABLE LAW, INCLUDING LIABILITY FOR DEATH OR PERSONAL INJURY.

Other damages

SUBJECT TO THE PARAGRAPH ABOVE, AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, Luca Bloom AND ITS OWNER, AGENTS, AND SUPPLIERS SHALL NOT BE LIABLE FOR ANY INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES, OR FOR ANY LOSS OF PROFITS, REVENUE, DATA, OR USE, ARISING OUT OF OR RELATING TO THIS PRODUCT OR THESE TERMS.

SUBJECT TO THE PARAGRAPH ABOVE, AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, Luca Bloom'S TOTAL AGGREGATE LIABILITY ARISING OUT OF OR RELATING TO THIS PRODUCT SHALL NOT EXCEED THE AMOUNT YOU ACTUALLY PAID FOR THE PRODUCT.

11. Severability & Entire Agreement

Version of these terms

These terms may be updated from time to time. The version that applies to you is the version in force on the date you purchased the product; its effective date appears at the end of this document. An updated version applies only to purchases made on or after its own effective date.

If any provision of these terms is held to be invalid, illegal, or unenforceable, that provision shall be modified to the minimum extent necessary to make it enforceable, or, if modification is not possible, severed. All remaining provisions shall continue in full force and effect.

These terms, together with the online manual referenced on the enclosed card, constitute the entire agreement between you and Luca Bloom regarding this product. No statement by any retailer, reseller, or third party modifies these terms. This paragraph does not limit any right you have under the Australian Consumer Law.

12. Contact

For questions, defect reports, warranty claims, or any formal notice under these terms, contact us at:

Luca Bloom

S-Building 3F, 2-1-19 Roppongi, Minato-ku, Tokyo, Japan

Telephone: +81 70-3842-6430

Telephone enquiries are answered in Japanese only, on weekdays from 10:00 to 17:00 Japan Standard Time. For enquiries in English, please use email or the platform messaging system; we answer in English by those channels.

Email: lucabloom65@gmail.com

You may also contact us through the messaging system of the platform where you purchased this product.

Seller: Luca Bloom

Components manufactured in China. Assembled, cleaned and inspected in Japan.

Effective Date: September 1, 2026